

Dividella

Pharma Packaging Technology

The Dividella *NeoTOP 904*



Range of formats:

min. 60 x 45 x 17
max. 150 x 150 x 32

Pack contents:

Ampoules, vials, cartridges,
syringes, tubes etc.

Range of volumes:

Vials: 2 - 20ml
Ampoules: 1 - 30ml
other sizes and syringes to order.

Output:

up to 900 objects/min., or 90 packs/min.

Packaging:

Printed or plain carton blanks from flat die cut blanks.
Partition walls from flat die cut blanks.

Application:

Depending on the outfitting of the machine, fully automatic forming and erecting of NeoTOP partitioned carton fluted trays, NeoTOP open-top trays or NeoTOP cartons including integrated partition from flat blanks; cylindrical objects are then inserted.

Method of operation:

The partitioned cartons, trays and packs are formed fully automatically from two supplied flat blanks. The partitions are formed and glued into the erected pack or onto the carrier blank. The products are then laid in by an automatic integrated inserter. An automatic product presence checker and leaflet inserter, barcode reading systems and a hot-melt unit for tamper-evident packaging complete the machine system.

Additional equipment:

Leaflets fed from sheets, barcode reading system for packs and leaflet, product and leaflet presence checking system, hot foil system, labelling on the inside lid of the pack, lasercoding, inserter for booklets.

Format change:

10 - 45 min (2 Persons) / no tool required

Machine size:

aprox. 9000 x 3300 x 2100 mm.
(depending on configuration)

